

# Conductivity standards



## Conductivity standards

Reagecon offer one of the largest ranges of certified conductivity standards with values ranging from 1.3  $\mu\text{S}/\text{cm}$  to 500,000  $\mu\text{S}/\text{cm}$ .

- Fully traceable – All values are directly traceable to NIST Certified Reference Material.
- Certified Accuracy – Range has an accuracy of  $\pm 1.0\%$  (except 1.3  $\mu\text{S}/\text{cm}$ )
- Guaranteed Stability
- Matrix Matched
- USP 24 compliant values.
- All Standards come complete with Certificate of Analysis detailing; test method used, result of test, traceability information, Lot No., Expiry date and comprehensive Measurement of Uncertainty Data.
- Temperature dependence data, Lot No. and expiry date are also printed on product labels.

Cat. No.	Code	Description	Pack Size
CO/61002	CSKC13	1.3 $\mu\text{S}/\text{cm}$	1 x 250mL
CO/61004	CSKC136	1.3 $\mu\text{S}/\text{cm}$	6 x 250mL
CO/61006	CSKC5	5 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61008	CSKC10	10 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61010	CSKC20	20 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61012	CSKC50	50 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61014	CSKC84	84 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61016	CSKC100	100 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61018	CSKCS	147 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61020	CSKC200	200 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61022	CSKC500	500 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61024	CSKC1000	1000 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61026	CSKCL	1413 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61028	CSKC10M	10,000 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61030	CSKC12880	12,880 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61032	CSKC20M	20,000 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61034	CSKC50M	50,000 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61036	CSKC100M	100,000 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61038	CSKC150M	150,000 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61040	CSKC200M	200,000 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61042	CSKC300M	300,000 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61044	CSKC350M	350,000 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61046	CSKC450M	450,000 $\mu\text{S}/\text{cm}$	1 x 500mL
CO/61048	CSKC500M	500,000 $\mu\text{S}/\text{cm}$	1 x 500mL



CO/610



CO/31120

## HANNA Calibration Conductivity/TDS/Salinity Solutions



Cat. No.	Code	Description
CO/31110	HI-7033M	84 $\mu\text{S}/\text{cm}$ Calibration solution, 230ml
CO/31112	HI-7033L	84 $\mu\text{S}/\text{cm}$ Calibration solution, 500ml
CO/31116	HI-7031M	1413 $\mu\text{S}/\text{cm}$ Calibration solution, 230ml
CO/31118	HI-7031L	1413 $\mu\text{S}/\text{cm}$ Calibration solution, 500ml
CO/31120	HI-70031P	1413 $\mu\text{S}/\text{cm}$ Calibration solution, 25 x 20ml sachets
CO/31122	HI-7030M	12.88 $\mu\text{S}/\text{cm}$ Calibration solution, 230ml
CO/31124	HI-7030L	12.88 $\mu\text{S}/\text{cm}$ Calibration solution, 500ml
CO/31126	HI-70030P	12.88 $\mu\text{S}/\text{cm}$ Calibration solution, 25 x 20ml sachets
CO/31128	HI-7032M	1382ppm Calibration solution, 230ml
CO/31130	HI-7032L	1382ppm Calibration solution, 500ml
CO/31132	HI-70032P	1382ppm Calibration solution, 25 x 20ml sachets
CO/31134	HI-7037L	100% NaCl Seawater standard solution, 500ml
CO/31136	HI-70442M	1500ppm Calibration solution, 230ml



CO/31134