

Indicator Papers

Whatman

Indicator paper, dispenser type

Supplied in 6mm x 5m reels with colour coded comparison chart printed on the dispenser.

Cat. No.	Code	Description
IP/17002	2600-100A	Full range pH 1.0 -14.0
IP/17004	2600-101A	Narrow range pH 0.5 - 5.5
IP/17006	2600-102A	Narrow range pH 4.0 - 7.0
IP/17008	2600-103A	Narrow range pH 6.0 - 8.0
IP/17010	2600-104A	Narrow range pH 8.0 -10.0
IP/17012	2600-201A	Litmus blue
IP/17014	2600-202A	Litmus red
IP/17016	2602-500A	Potassium Iodide



IP/170

Indicator papers

Supplied in boxes containing 10 books, each book contains 20 leaves.

Cat. No.	Description
IP/97002	Wide range universal pH 1-14
IP/97004	Wide range universal pH 1-11
IP/97006	Narrow range pH 1.0-3.5
IP/97008	Narrow range pH 3.6-5.1
IP/97010	Narrow range pH 5.2-6.7
IP/97012	Narrow range pH 6.8-8.3
IP/97014	Narrow range pH 8.4-10.0
IP/97016	Litmus red
IP/97018	Litmus blue
IP/97020	Litmus neutral
IP/97022	Cobalt chloride
IP/97024	Lead acetate
IP/97026	Starch
IP/97028	Starch/Iodide
IP/97030	Turmeric (20 strips)



IP/97004

Indicator test strips, pH

Made from non bleed test paper mounted on plastic wand. Includes comparison chart giving readings as detailed below. Supplied in packs of 50 strips.

Cat. No.	Range	Chart readings
IP/97102	0 to 1.5	pH 0, 0.3, 0.5, 0.8, 1.0 & 1.5
IP/97104	0 to 2.5	pH 0, 0.5, 1.0, 1.5, 2.0 & 2.5
IP/97106	3 to 6	pH 3.0, 3.5, 4, 4.5, 5, 5.5 & 6.0
IP/97108	5 to 9	pH 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, & 9.0
IP/97110	7 to 10	pH 7.0, 7.5, 8, 8.5, 9, 9.5 & 10.0
IP/97112	10 to 13	pH 10.0, 10.5, 11.0, 11.5, 12.0, 12.5 & 13.0
IP/97114	11.5 to 14	pH 11.5, 12.0, 12.5, 13.0, 13.5 & 14.0
IP/97116	0 to 14	pH 0,1,2,3,4,5,6,7,8,9,10,11,12,13 & 14



IP/9710 series

Indicator papers, Chlorine test strips

Complete with comparison chart giving readings as detailed below. Supplied in packs of 100 strips.

Cat. No.	Test paper	Chart readings
IP/97120	Low-level chlorine	0, 1.0, 3.0, 5.0 and 10 mg/L (ppm)
IP/97122	Chlorine	0, 25, 50, 100 and 300 mg/L (ppm)



IP/9712

Indicator papers, water hardness strips

Used to monitor water quality in laundry operations as well as home and laboratory use. Includes comparison chart giving readings as detailed below. Supplied in packs of 50 strips.

Cat. No.	Chart readings
IP/97202	0, 5, 125, 250 and 500 mg/L (ppm).



IP/97202

Indicator papers, general purpose strips

Used in many industrial as well as laboratory applications. Includes comparison chart giving readings as detailed below. Supplied in packs of 50 strips.

Cat. No.	Description	Range	Chart readings
IP/97302	Sulphite test	10-500mg/L	10, 50, 100, 250 & 500mg/L (ppm)
IP/97304	Nitrite test	0-25mg/L	0,0.5, 1, 5,10 & 25mg/L (ppm)
IP/97306	Nitrate test	0-500mg/L	0,10,25,50,100, 250 & 500mg/L
IP/97310	Peroxide test	1-100mg/L	0, 1, 3, 10, 30 & 100mg/L (ppm)



IP/97306