



SR/760



SR/761



SR/762

## Stoppers

### Stoppers, natural red rubber

Available as either solid stoppers, stoppers with one hole or stoppers with two holes.

Standard holes are 3.0mm diameter in the 5 to 8mm base diameter single hole stoppers and 5.0mm diameter in all other sizes of one and 2 hole stoppers.

Cat. No. Solid	Cat. No. One hole	Cat. No. Two holes	Bottom dia. mm	Top dia. mm	Length mm	Pk Qty
SR/76002	~	~	3	4.5	12	50
SR/76004	~	~	4	5.8	14	50
SR/76006	SR/76106	~	5	7	16	50
SR/76008	SR/76108	~	6	8	16	50
SR/76010	SR/76110	~	7	9	16	20
SR/76012	SR/76112	~	8	10.5	20	20
SR/76014	SR/76114	~	9	11.5	20	20
SR/76016	SR/76116	~	10	12.5	20	20
SR/76018	SR/76118	~	11	14	24	20
SR/76020	SR/76120	~	13	16	24	20
SR/76022	SR/76122	SR/76222	15	18	24	20
SR/76024	SR/76124	SR/76224	17	20.3	26	20
SR/76026	SR/76126	SR/76226	19	22.5	28	20
SR/76028	SR/76128	SR/76228	21	24.5	28	20
SR/76030	SR/76130	SR/76230	23	26.5	28	20
SR/76032	SR/76132	SR/76232	25	28.5	28	20
SR/76034	SR/76134	SR/76234	27	31	32	10
SR/76036	SR/76136	SR/76236	29	33	32	10
SR/76038	SR/76138	SR/76238	31	36	35	10
SR/76040	SR/76140	SR/76240	33	38.5	38.5	10
SR/76042	SR/76142	SR/76242	35	40.5	38.5	10
SR/76044	SR/76144	SR/76244	37	42.5	38.5	10
SR/76046	SR/76146	SR/76246	39	44.5	38.5	10
SR/76048	SR/76148	SR/76248	41	47	42	5
SR/76050	SR/76150	SR/76250	43	49	42	5
SR/76052	SR/76152	SR/76252	45	51	42	5
SR/76054	SR/76154	SR/76254	49	55	42	5
SR/76056	SR/76156	SR/76256	53	59.5	45.5	2
SR/76058	SR/76158	SR/76258	57	64	49	2
SR/76060	SR/76160	SR/76260	62	69	49	2
SR/76062	SR/76162	SR/76262	67	74	49	2
SR/76064	SR/76164	SR/76264	71	78	49	1
SR/76066	SR/76166	SR/76266	75	85	54	1
SR/76068	SR/76168	SR/76268	79	89	54	1
SR/76070	SR/76170	SR/76270	86	96	54	1
SR/76072	SR/76172	SR/76272	92	102	54	1
SR/76074	SR/76174	SR/76274	99	109	54	1
SR/76076	SR/76176	SR/76276	105	116.5	63.5	1

Stoppers are also available in grey Neoprene - please ask for details.

Assorted pack of 67 solid or natural rubber stoppers, sizes included, 5 of each; 7, 8, 9, 10, 11, 13, 15, 17, 19, 21, 23 and 2 of each 25, 27, 29, 31, 33, 35mm

Cat. No.	Description
SR/76080	Assorted stoppers, red rubber, solid, pack of 67
SR/76082	Assorted stoppers, red rubber, 1 hole, pack of 67
SR/76084	Assorted stoppers, red rubber, 2 holes, pack of 48 (5 of each sizes 15-29 and 2 each 31-37)