



SY/920

Hypodermic, syringe, all glass, Luer fitting, without needle

Cat. No.	Capacity, ml
SY/92002	1
SY/92004	2
SY/92006	5
SY/92008	10
SY/92010	20



SY/930

Glass, syringe, metal Luer lock fitting, without needle

Cat. No.	Capacity, ml
SY/93012	1
SY/93014	2
SY/93016	5
SY/93018	10



SY/931

Glass, syringe, Luer fitting with metal tip, without needle

Cat. No.	Capacity, ml
SY/93104	2
SY/93106	5
SY/93108	10
SY/93110	20
SY/93112	50
SY/93114	100

Gastight syringes, 1700 series

- Capacities 10-500µl
- Cemented needles (LTN), removable needles (RN), luer tip (LT)
- Precision-machined PTFE-tipped plungers
- Replaceable plungers



Syringe range

Cat. No.	Description
Luer tip	
SY/90102	Syringe, gas tight 1701 luer tip, 10µl
SY/90104	Syringe, gas tight 1702 luer tip, 25µl
SY/90106	Syringe, gas tight 1705 luer tip, 50µl
SY/90108	Syringe, gas tight 1710 luer tip, 100µl
SY/90110	Syringe, gas tight 1725 luer tip, 250µl
SY/90112	Syringe, gas tight 1750 luer tip, 500µl
Cemented needle	
SY/90202	Syringe, gas tight 1701 cemented needle, 10µl
SY/90204	Syringe, gas tight 1702 cemented needle, 25µl
SY/90206	Syringe, gas tight 1705 cemented needle, 50µl
SY/90208	Syringe, gas tight 1710 cemented needle, 100µl
SY/90210	Syringe, gas tight 1725 cemented needle, 250µl
Removable needle	
SY/90302	Syringe, gas tight 1701 removable needle, 10µl
SY/90304	Syringe, gas tight 1702 removable needle, 25µl
SY/90306	Syringe, gas tight 1705 removable needle, 50µl
SY/90308	Syringe, gas tight 1710 removable needle, 100µl
SY/90310	Syringe, gas tight 1725 removable needle, 250µl
SY/90312	Syringe, gas tight 1750 removable needle, 500µl

Gastight syringes, 1000 series

- Capacities 1 - 100ml
- Cemented needles (LTN), removable needles (RN), luer tip (LT)
- Precision-machined PTFE-tipped plungers
- Replaceable plungers

Cat. No.	Description
Luer tip	
SY/90402	Syringe, gas tight 1001 luer tip, 1ml
SY/90404	Syringe, gas tight 1002 luer tip, 2.5ml
SY/90406	Syringe, gas tight 1005 luer tip, 5ml
SY/90408	Syringe, gas tight 1010 luer tip, 10ml
Cemented needle	
SY/90502	Syringe, gas tight 1001 cemented needle, 1ml
SY/90504	Syringe, gas tight 1002 cemented needle, 2.5ml
SY/90506	Syringe, gas tight 1005 cemented needle, 5ml
SY/90508	Syringe, gas tight 1010 cemented needle, 10ml
Removable needle	
SY/90602	Syringe, gas tight 1001 removable needle, 1ml
SY/90604	Syringe, gas tight 1002 removable needle, 2.5ml
SY/90606	Syringe, gas tight 1005 removable needle, 5ml
SY/90608	Syringe, gas tight 1010 removable needle, 10ml